

**EMS Order form**  
**Multi-Surface Flexible Inner Socket**

**Send fiberglass cast & order form to:**

**NETTwork EMS**  
**40091 County Road 1**  
**Rice, MN 56367**

Date Ordered: \_\_\_\_\_

Purchase Order Number: \_\_\_\_\_

Ordered by: \_\_\_\_\_

Bill to Address: \_\_\_\_\_

Prosthetist: \_\_\_\_\_

Phone: \_\_\_\_\_

Email Address: \_\_\_\_\_ Billing Email Address: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

**Shipping preference:** UPS ground (no charge) \_\_\_\_\_ Other (please specify) \_\_\_\_\_

***Patient Information:***

Last Name: \_\_\_\_\_ First Name: \_\_\_\_\_

Patient K-Level: K1: \_\_\_\_\_ K2: \_\_\_\_\_ K3: \_\_\_\_\_ K4: \_\_\_\_\_

Left: \_\_\_\_\_ Right: \_\_\_\_\_

Transtibial: \_\_\_\_\_ Transfemoral: \_\_\_\_\_

***Items to be Ordered Information: Pricing below valid when ordering directly from NETTwork EMS***

EMS "Multi-Surfaced Flexible Inner Socket": \_\_\_\_\_ (Your Cost \$950.00)

Corrugation Reductions: \_\_\_\_\_ (mark areas within the cast radially including distal end as an option)

Standard Test-Socket with Distal end Air Chamber: \_\_\_\_\_ (Your Cost \$225.00 when ordered with EMS)

Test-Socket without Distal end Air Chamber: \_\_\_\_\_ (Your Cost \$225.00 when ordered with EMS)

Additional EMS "Multi-Surfaced Socket": \_\_\_\_\_ (Your Cost \$475.00 when ordered with the first EMS)

***Prosthesis Application Information:***

Elevated Vacuum: \_\_\_\_\_ Suction type System: \_\_\_\_\_ Pin Type system: \_\_\_\_\_

To re-order EMS without changes, indicate the Inner Sockets Serial Number: \_\_\_\_\_

***Payment Information:***

Full payment due Net 30. Now accepting credit cards! (subject to 3.5% processing fee)

**Warranty information:**

The EMS Multi-Surfaced Inner Flexible Socket is the most reliable way to ensure a patient's prosthetic limb remains securely in place. While our proprietary material and manufacturing process has proven to produce durable inner sockets second to none, we understand that sometimes things happen. Every EMS inner socket comes with a six (6) month warranty which begins the day your order ships.

If you (patient or prosthetist) suspect a material failure or fabrication defect, call us immediately at 866-654-0858

You may be asked to provide the following: purchase order number, serial number, date suspected defect was observed, and photos. Once the warranty exchange is approved, you will be shipped a replacement EMS socket along with a return shipping label to send us the defective socket.

We value each and every customer & we appreciate the opportunity to assist in providing a better quality of life for amputees.

*Prosthetist's Notes to Manufacturer:*

---

---

---

---

For additional information, please visit us at:

[www.NETTworkMFG.com/EMS](http://www.NETTworkMFG.com/EMS)